



Recombinant Human cytomegalovirus G-protein coupled receptor homolog US27 (US27)

Product Code	CSB-CF357733HWV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09703
Product Type	Transmembrane Protein
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Sequence	MTTSTNNQTLTQVSNMTNHTLNSTEIYQLFEYTRLGVWLMCIVGTFLNVLVITTI LYYRR KKKSPSDTYICNLAVADLLIVVGLPFFLEYAKHHPKLSREVVCSGLNACFYICLF AGVCF LINLSMDRYCVIVWGVELNRVRNNKRATCWWVIFWILAVLMGMPHYLMYSHTN NECVGEF ANETSGWFVPVFLNTKVNICGYLAPIALMAYTYNRMVRFIINYVGKWHMQLHLVL LVVVVS FASFWFPFNALFLESIRLLAGVYNDTLQNVIIIFCLYVGQFLAYVRACLNPGIYIL VGTQ MRKDMWTTLRVFACCCVKQEIPYQDIDIELQKDIQRRAKHTKRTHYDRKNAPM ESGEEEF LL
Source	in vitro E.coli expression system
Target Names	US27
Protein Names	Recommended name: G-protein coupled receptor homolog US27 Alternative name(s): HHRF2
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.