



Recombinant Human Olfactory receptor 5AU1 (OR5AU1)

Product Code	CSB-CF839884HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGC0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MTEFHLQSQMPSIRLIFRRLSLGRIKPSQSPRCSTSFMVVPSFSIAEHWRRMK GANLSQG MEFELLGLTTDPQLQRLLFVVFLGMYTATLLGNLVMFLLIHVSATLHTPMYSLLK SLSFL DFCYSSTVVPQTLVNFLAKRKVISYFGCMTQMFFYAGFATSECYLIAAMAYDR YAAICNP LLYSTIMSPEVCASLIVGSYSAGFLNSLIHTGCIFSLKFCGAHVVTHFFCDGPPIL SLSC VDTSLCEILLFIFAGFNLLSCTLTILISYFLILNTILKMSSAQGRFKAFSTCASHLTA IC LFFGTTLFMYLRPRSSYSLTQDRTVAVIYTVVIPVLNPLMYSLRNKDVKKALIKV WGRKT ME
Source	in vitro E.coli expression system
Target Names	OR5AU1
Protein Names	Recommended name: Olfactory receptor 5AU1 Alternative name(s): Olfactory receptor OR14-38
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.