



Recombinant Pig G-protein coupled receptor 4 (GPR4)

Product Code	CSB-CF009813PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50132
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MGNGTWEGCHVDSRVDHLFPPSLYIFVIGVGLPTNCLRLWAAAYRQVRQRNEL GVYLMNLS IADLLYICTLPLWVDYFLHHDNWIHGPGSCKLFGFIFYTNIYISIAFLCCISVDRYL AVA HPLRFARLRRVKTAVAVSSVWATELGANSVPLFHDELFRDRYNHTFCFEKFP MEGWVAW MNLYRVFVGFLFPWALMLLSYRGILRAVRGSVSTERQEKAKIKRLALSIAIVLV CFAPY HVLLLSRSAVYLGHPWDCGFEERVFSAYHSSLAFTSLNVCVADPILYCLVNEGA RSDVAKA LHNLLRFLTSDKPKQEMASASLTLDTPLTSKRNSMARAVAAAGWHLRPPRDQV QLKMLPPP AP
Source	in vitro E.coli expression system
Target Names	GPR4
Protein Names	Recommended name: G-protein coupled receptor 4 Alternative name(s): G-protein coupled receptor 19
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.