



Recombinant Serpentine receptor class beta-18 (srb-18)

Product Code	CSB-CF519421CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O16444
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	<p>MTSEEIHAAGLCSVVPTIFQNCNDWISKNGTCYQDSCCMSLETATSWYYRLA QFSHVVF SFMGLIIVVVYILRYRSRHILPENVRVLVDFMLLFIVAHSIDMIVLHIYHIIQSFQANI S DPCFVREKVSFCAPFRYTFSFCSMGLAICTYCIYIDRLACAYYKNYTKHQRLILA AQICQ LIVISSLIIVVYRNEEPNTYLLSCLNVPVASVEDMAKATIAVFPINFICFFLSIGLF RH FKKKEEGSRFDIVRHFTASVDVESSEFLFRTTGTQAALMALFSVASLLMRLVYN FLPRQV GLTIATLSYIMSIYCFTVPLVIVKCVQKTSALRKSRISSHVGLKAMGVEGASNYF EMMKS QWE</p>
Source	in vitro E.coli expression system
Target Names	srb-18
Protein Names	Recommended name: Serpentine receptor class beta-18 Short name= Protein srb-18
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.