



Recombinant Oncorhynchus mykiss C3a anaphylatoxin chemotactic receptor (c3ar1)

Product Code	CSB-CF650496OEI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2WED0
Product Type	Transmembrane Protein
Immunogen Species	Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri)
Sequence	MGDNMDFSEHYGNFSENYVTESYGEFDLYDPLNETSLSEQGHRSIWVLSIVL CSIACVL GITGNAFVIWIAGVKMKRTVNTIWFVNLAADLLCCVSIPFSIADIILNSHWPYGE AMCK ILPSMVVLNMFASVFTLVLI SLDRFALVILPVWAQNHR SITLAWLLCGLVWVLGL LLSLP SMIYREIVVHDDMNITLCIYNHLQDKTEGNQSAIKAIHVTRLILGFLIPLLVIACYL LI GRRVSSGRFKSQRAFQIILVVVTTFFVCWLPYHVIGLVIEYGKEASQVMARALD PLAISL AYNNSCLNPVLYVFMGQDFKERVVRVSLRKIFEKVFSEDVTLRSSVYSKGQSQL SRATNSS EAQV
Source	in vitro E.coli expression system
Target Names	c3ar1
Protein Names	Recommended name: C3a anaphylatoxin chemotactic receptor Short name= C3AR Short name= C3a-R
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.