



Recombinant Bovine Long-wave-sensitive opsin 1 (OPN1LW)

Product Code	CSB-CF858376BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BGI7
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	<p>MAHAWGPQRLAGGQPQANFEESTQGSIFTYTNSNSTRDPFEGPNYHIAPRW VYHLTSAWM VFVVIASVFTNGLVLAATMRFKCLRHPLNWILVNLAIDLAEETIIASTISVVNQMY GYFV LGHPLCVVEGYTVSLCGITGLWLSLAIISWERWMVVKPFNGVRFDAKLAITGIA FSWIWA AVWTAPPPIFGWSRYWPHGLKTSCGPDVFSGSSYPGVQSYMIVLMITCCFIPLS VIILCYL QVWLAIRAVAKQQKESESTQKAEKEVTRMVMVMIFAYCLCWGPYTTFFACFAA AHPGYAFH PLVAALPAYFAKSATIYNPIIYVFMNRQFRNCILQLFGKKVDDSSSELSSVSKTEA SSVSS VSPA</p>
Source	in vitro E.coli expression system
Target Names	OPN1LW
Protein Names	Recommended name: Long-wave-sensitive opsin 1 Alternative name(s): Red cone photoreceptor pigment Red-sensitive opsin
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.