



Recombinant Mouse Growth hormone secretagogue receptor type 1 (Ghsr)

Product Code	CSB-CF858781MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q99P50
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MWNATPSEEPEPNVTLDLDWDASPGNDSLDELLPLFPAPLLAGVTATCVALF VVGISGN LLTMLVVSRLFRELRTTTNLYLSSMAFSDLLIFLCMPLDLVRLWQYRPWNFGDLL CKLFQF VSESCYATVLTITALSVERYFAICFPLRAKVVVTKGRVKLVILVIWAVAFCSAG PIFVL VGVEHENGTDPRDTNECRATEFAVRSGLLTVMVWVSSVFFFLPVFCLTVLYSL IGRKLWR RRGDAAVGSSLRDQNHKQTVKMLAVVVFAFILCWLPFHVGRYLFKSKSFEPGS LEIAQISQ YCNLVSFVLFYLSAAINPILYNIMSKKYRVAVFKLLGFESFSQRKLSTLKDESSR AWTKS SINT
Source	in vitro E.coli expression system
Target Names	Ghsr
Protein Names	Recommended name: Growth hormone secretagogue receptor type 1 Short name= GHS-R Alternative name(s): GH-releasing peptide receptor Short name= GHRP Ghrelin receptor
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.