



Recombinant Human Somatostatin receptor type 5 (SSTR5)

Product Code	CSB-CF022728HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35346
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEPLFPASTPSWNASSPGAASGGGDNRTLVGPPSAGARAVLVPVLYLLVCA AGLGGNTL VIYVVLRFKMKTVTNIYILNLAVADVLYMLGLPFLATQNAASFVWPGPVLCLRLV MTLDG VNQFTSVFCLTVMSVDRYLAVVHPLSSARWRRPRVAKLASAAAWVLSLCMSL PLLVFADV QEGGTCNASWPEPVGLWGAVFIIYTAVLGGFFAPLLVICLCYLLIVVKVRAAGVR VGCVRR RSERKVTRMVLVVVLVFAGCWLPFFFTVNIVNLAVALPQEPASAGLYFFVILSY ANSCAN PVLYGFLSDNFRQSFQKVLCLRKGSAGKADADATEPRPDRIRQQQEATPPAHR AAANGLMQ TSKL
Source	in vitro E.coli expression system
Target Names	SSTR5
Protein Names	Recommended name: Somatostatin receptor type 5 Short name= SS-5-R Short name= SS5-R Short name= SS5R
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.