



Recombinant Rat Lysophosphatidic acid receptor 1 (Lpar1)

Product Code	CSB-CF013047RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61794
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MAAASTSSPVISQPQFTAMNEQQCFYNESIAFFYNRSGKYLATEWNTVSKLVM GLGITVC VFIMLANLLVMVAIYVNRFFHFPIYYLMANLAAADFFAGLAYFYLMFNTGPNTR RLTVST WLLRQGLIDTSLTASVANLLAIAIERHITVFRMQLHTRMSNRRVVVVIVVIWTMAI VMGA IPSVGWNCICDIDHCSNMAPLYSDSYLVFWAIFNLVTFVVMVLYAHIFGYVRQ RTMRMS RHSSGPRRNRDTMMSLLKTVVIVLGAFIVCWTPGLVLLLLDVCCPQCDVLAYE KFFLLLA EFNSAMNPIIYSYRDKEMSATFRQILCCQRNENPNGPTEGSDRSASSLNHTILA GVHSND HSVV
Source	in vitro E.coli expression system
Target Names	Lpar1
Protein Names	Recommended name: Lysophosphatidic acid receptor 1 Short name= LPA receptor 1 Short name= LPA-1 Alternative name(s): Lysophosphatidic acid receptor Edg-2
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.