



Recombinant Rabbit Medium-wave-sensitive opsin 1 (OPN1MW)

| | |
|--------------------------|---|
| Product Code | CSB-CF016352RB |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O18910 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Oryctolagus cuniculus (Rabbit) |
| Sequence | <p>MTQPWGPQMLAGGQPPEshedstQASIFTYtNSNSTRGPFEGPNFHIAPRW VYHLTSAWM ILVVIASVFTNGLVLVATMRFKCLRHPLNWILVNLAVADLAETVIASISVVNQFY GYFV LGHPLCVVEGYTVSLCGITGLWSLAIISWERWLVVCKPFGNVRFDKLAIA GIA FSWIWA AVWTAPPIFGWSRYWPYGLKTSCGPDVFSGTSYPGVQSYMMVLMVTCCIIPL SVIVLCYL QVWMAIRTVAKQKQKESSTQKAEKEVTRMVVVMVFAYCLCWGPYTFACFAT AHPGYSFH PLVAAIPSYFAKSATIYNPIIYVFMNRQFRNCILQLFGKKVEDSSELSSASRTEAS SVSS VSPA</p> |
| Source | in vitro E.coli expression system |
| Target Names | OPN1MW |
| Protein Names | Recommended name: Medium-wave-sensitive opsin 1 Alternative name(s): Green cone photoreceptor pigment Green-sensitive opsin Medium wavelength-sensitive cone opsin |
| Expression Region | 1-364 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |