



Recombinant Mouse Probable G-protein coupled receptor 142 (Gpr142)

Product Code	CSB-CF742517MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TQN9
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MHLNSNPNSYICDAYQHADLLWSLSPHVLTKAVQPQVTLPTVNGSNPRYDG VDGHWPEP PERSPCVAGIIPVIYYSVLLSLGLPVALARLAARTRKPSYHYLLALTASDIVTQVII VFV GFLQGAFLARQVPQAVVRTANILEFAANHASVWIAVLFTVDRYNALCRPLRH RATSSPG RTHRAIAAVIGVTLTGIPFYWWLDVWRDADPPSTMDKLLKWAHCLIVYFIPCN VFLVTN SAIILRLRKRGRGLRPLVSKSTAILLGVTSLFALLWAPRIIVMLYHLYVAPVHRD WRVH LALDIANMLAMLNTEVNFGLYCFISKTFRATVRQVICDVHMACALKSQPKQTVV ELMLKS VGTEL
Source	in vitro E.coli expression system
Target Names	Gpr142
Protein Names	Recommended name: Probable G-protein coupled receptor 142
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.