



Recombinant Human 5-hydroxytryptamine receptor 1E (HTR1E)

Product Code	CSB-CF010885HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28566
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNITNCTTEASMAIRPKTITEKMLICMTLVVITTLTLLNLAVIMAIGTTKKLHQPA NYL ICSLAVTDLLVAVLVMPLSIYIVMDRWKLG YFLCEVWLSVDMTCCTCSILHLCVI ALDR YWAITNAIEYARKRTAKRAALMILTVWTISIFISMPPLFWRSRRLSPPPSQCTI QHDHV IYTIYSTLGAFYIPLTLILILYYRIYHAAKSLYQKRGSSRHLSNRSTDSQNSFASCK LTQ TFCVSDFFSTSDPTTEFEKFKHASIRIPFDNDLDHPGERQQISSTRERKAARILGL ILGAF ILSWLPFFIKELIVGLSIYTVSSEVADFLTWLGYVNSLINPLLYTSFNEDFKLAFKK LIR CREHT
Source	in vitro E.coli expression system
Target Names	HTR1E
Protein Names	Recommended name: 5-hydroxytryptamine receptor 1E Short name= 5-HT-1E Short name= 5-HT1E Alternative name(s): S31 Serotonin receptor 1E
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.