



Recombinant Human P2Y purinoceptor 4 (P2RY4)

Product Code	CSB-CF017335HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51582
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MASTESSLLRSLGLSPGPGSSEVELDCWFDEDFKFILLPVSYAVVFLGLGLNA PTLWLF IFRLRPWDATATYMFHLALSDTLYVLSLPTLIYYAAHNHWPFGEICKFVRFLF YWNLY CSVLFLTCISVHRYLGICHPLRALRWGRPRLAGLLCLAVWLVVAGCLVPNLFFV TTSNKG TTVLCHDTRPEEFDHYVHFSSAVMGLLFGVPCLVTLVCYGLMARRLYQPLPG SAQSSSR LRSLRTIAVVLTVFAVCFVPFHITRTIYYLARLLEADCRVLNIVNVVYKVTRPLAS ANSC LDPVLYLLTGDKYRRQLRQLCGGGKQPRTAASSLALVSLPEDSSCRWAATP QDSSCSTP RADRL
Source	in vitro E.coli expression system
Target Names	P2RY4
Protein Names	Recommended name: P2Y purinoceptor 4 Short name= P2Y4 Alternative name(s): P2P Uridine nucleotide receptor Short name= UNR
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.