



Recombinant Bovine C-X-C chemokine receptor type 3 (CXCR3)

Product Code	CSB-CF683649BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5MD61
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MVPEMSERQEFQASDFAYLLENSSYDYGENETYFCCTSPPCPQDFSLNFDRT FLPVLYSL LFVLGLLGNGIVAVVLLSQRAALSSTDTFLLHLAVADALLVLTPLWAVDAAIQW VFGSG LCKVAGALFNINFYAGALLLACISFDRLSIVHATQLYRRGPPTRVALTCVAVW GLCLLF ALPDFIFLSSHHDNRLNATHCQYNFPQEGHTALRILQLVAGFLLPLLVMAYCYA RILAVL LVSQRRLRAMRLVVVVVAFALCWTPYHLVVLVDTLMDLGALARNCGRES SVDIAKSV TSGMGYMHCCLNPLLYAFVGVKFRERMWVLLVRLGCPDQRCHQRQPSASRR ESSWSETTE ASYSGL
Source	in vitro E.coli expression system
Target Names	CXCR3
Protein Names	Recommended name: C-X-C chemokine receptor type 3 Short name= CXC-R3 Short name= CXCR-3 Alternative name(s): Interferon-inducible protein 10 receptor Short name= IP-10 receptor CD_antigen= CD183
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.