



Recombinant Pig Growth hormone secretagogue receptor type 1 (GHSR)

Product Code	CSB-CF853264PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95254
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MWNATPSEEPGPNLTLPDLGWDAPPENDSLVEELLPLFPTPLLAGVTATCVAL FVVGIAQ NLLTMLVVSFRFREMRTTNTNLYLSSMAFSDLLIFLCMPLDLFRLWQYRPWNLGN LLCKLFQ FVSESCYATVLTITALSVERYFAICFPLRAKVVVTKGRVKLVILVIWAVAFCSAG PIFV LVGVEHDNGTDPRDTNECRATEFAVRSGLLTVMVWVSSVFFFLPVFCLTVLYS LIGRKLW RRKRGEAAVGSSLRDQNHKQTVKMLAVVVFILCWLPFHVGRYLFSKSLEP GSVEIAQI SQYCNLVSVFLFYLSAAINPILYNIMSKKYRVAVFKLLGFEPFSQRKLSTLKDES SRAWT ESSINT
Source	in vitro E.coli expression system
Target Names	GHSR
Protein Names	Recommended name: Growth hormone secretagogue receptor type 1 Short name= GHS-R Alternative name(s): GH-releasing peptide receptor Short name= GHRP Ghrelin receptor
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.