



Recombinant Guinea pig 5-hydroxytryptamine receptor 1F (HTR1F)

Product Code	CSB-CF010886GU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O08890
Product Type	Transmembrane Protein
Immunogen Species	Cavia porcellus (Guinea pig)
Sequence	MDFLNSSDQNLTSSELLHRMPSKILVSLTSLGLALMTTINSLVIAAIIIVTRKLHH PANY LICSLAVTDFLVAVLVMPFSIVYIVRESWIMGQVLCDIWLSVDIICCTCSILHLSAI ALD RYRAITDAVEYARKRTPKQAGIMITIVWIISVFISMPPLFWRHQGTSRDDECIKH DHIV STIYSTFGAFYIPLVLILILYYKIYKAAKTLYHKRQASRIAKEELNGQVLLESGEKSI KM VSTTYVPEKSLSDPSTDFDKIHNTVKSPRCKLRHEKSWRRQKISGTRERKAAT TLGLILG AFVICWLPFFVKELVVNVCEKCKISEEMANFLAWLGYLNSLINPLIYTIFNEDFKK AFQK LVRCQY
Source	in vitro E.coli expression system
Target Names	HTR1F
Protein Names	Recommended name: 5-hydroxytryptamine receptor 1F Short name= 5-HT-1F Short name= 5-HT1F Alternative name(s): Serotonin receptor 1F
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.