



Recombinant Rat C-X-C chemokine receptor type 3 (Cxcr3)

Product Code	CSB-CF884915RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JII9
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MYLEVSERQVLDASDIAFLLENSTSPYDYGENESDFSDSPPCPQDFSLNFDRT FLPVLYS LLFLLGLLGNGAVAAVLLSQRTALSSTDTFLLHLAVADVLLVLTPLWAVDAAA QWVFGS GLCKVAGALFNINFYAGAFLLACISFDRLSIVHATQIYRRDPWVRVALTCIVVW GLCVL FALPDFIFLSASHDQRLNATHCQYNFPQVGRTRALRVLQLVAGFLMPLLVMAYC YAHILAV LLVSRGQRRFRAMRLVVVVVAVAVCWTPYHLVVLVDILMDVGVLARNGRE SHVDVAKS VTSGMGYMHCCLNPLLYAFVGVKFKEQMWMLLMRLGRSDQRGPQRQPSSS RRESSWSETT EASYLGL
Source	in vitro E.coli expression system
Target Names	Cxcr3
Protein Names	Recommended name: C-X-C chemokine receptor type 3 Short name= CXC-R3 Short name= CXCR-3 Alternative name(s): Interferon-inducible protein 10 receptor Short name= IP-10 receptor CD_antigen= CD183
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.