



Recombinant Bovine Somatostatin receptor type 2 (SSTR2)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-CF022725BO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P34993 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Sequence | MDLVSELNETQPWLTPFDLNGSVGAANISNQTEPYD LASNVVLTFIYFVVCII GLCGN TLVIYVILRYAKMKTITNIYILNLIAIADEL FMLGLPFLAMQVALVHWPFKGAICRVV MTV DGINQFTSIFCLTVMSIDRYLAVVHPIKSAKWRRPRTAKMINVAVWGV SLLVILPI MIYA GLRSNQWGRSSCTINWPGESGAWYTGFIYAFILGFLVPLTIICLCYLFIIIKVKSS GIR VGSSKRKKSEKKVTRMVSIVVAVFIFCWLPFYIFNVSSVSAISPTPALKGMFDF VVVLT YANSCANPILYAFLSDNFKKSFQNVLCLVKVSGTDDGERSDSKQDKSRLNETT ETQRTLL NGDLQTSI |
| Source | in vitro E.coli expression system |
| Target Names | SSTR2 |
| Protein Names | Recommended name: Somatostatin receptor type 2 Short name= SS-2-R Short name= SS2-R Short name= SS2R Alternative name(s): SRIF-1 |
| Expression Region | 1-368 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |