



Recombinant Dog Somatostatin receptor type 2 (SSTR2)

Product Code	CSB-CF676611DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q49LX6
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	MDMEYELLNESHTWVSPFDLDGSVVAANSSNQTEPYDLSNAVLTFIYFVV CIIGLCG NTLVIIYVILRYAKMKTITNIYILNLAIAGELFMLGLPFLAMQVALVHWPFGKAICRV VMT VDGINQFTSIFCLTVMSIDRYLAVVHPIKSAKWRRPRTAKMVNVAVWGVSLLVV LPIMIY AGLRSNQWGRSSCTINWPGESGTWYTGFIITYFILGFLVPLTIICLCYLFIIKVKKS SGI RVGSSKRKKSEKKVTRMVSIVVAVFIFCWLPFYIFNVSSVSVAISPTPALKGMF DLVVVL TYANSCANPILYAFLSDNFKKSFQNVLCLLVKVS GPDDGERSDSKQDKSRLNET TETQRTL LNGDLQTSI
Source	in vitro E.coli expression system
Target Names	SSTR2
Protein Names	Recommended name: Somatostatin receptor type 2 Short name= SS-2-R Short name= SS2-R Short name= SS2R
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.