



Recombinant Mouse Somatostatin receptor type 2 (Sstr2)

Product Code	CSB-CF022725MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30875
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEMSSEQLNGSQVWVSSPFDLNGSLGPSNGSNQTEPYDMTSNAVLTFIYFV VCVVGLCG NTLVIYVILRYAKMKTITNIYILNLAIADELFMLGLPFLAMQVALVHWPFGKAICRV VMT VDGINQFTSIFCLTVMSIDRYLAVVHPIKSAKWRRPRTAKMINVAVWCVSLLVIL PIMIY AGLRSNQWGRSSCTINWPGESGAWYTGFIYAFILGFLVPLTIICLCYLFIIKVKKS SGI RVGSSKRKKSEKKVTRMVSIVVAVFIFCWLPFYIFNVSSVSVAISPTPALKGMF DFVVIL TYANSCANPILYAFLSDNFKKSFQNVLCLLVKVSGTEDGERSDSKQDKSRLNET TETQRTL LNGDLQTSI
Source	in vitro E.coli expression system
Target Names	Sstr2
Protein Names	Recommended name: Somatostatin receptor type 2 Short name= SS-2-R Short name= SS2-R Short name= SS2R Alternative name(s): SRIF-1
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.