



Recombinant Human Prolactin-releasing peptide receptor (PRLHR)

Product Code	CSB-CF018726HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49683
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MASSTTRGPRVSDLFSGLPVAVTTPANQSAEASAGNGSVAGADAPAVTPFQS LQLVHQLK GLIVLLYSVVVVVGLVGNCLLVLIARVRRLLHNVTNFLIGNLALSDVLMCTACVP LTLAY AFEPRGWVFGGGLCHLVFFLQPVTVYVSVFTLTIIAVDRYVVLVHPLRRRISLR LSAYAV LAIWALSAVLALPAAVHTYHVELKPHDVRRLCEEFWGSQERQRQLYAWGLLLVT YLLPLL ILLSYVRVSVKLRNRVVPVCVTQSQADWDRARRRRTFCLLVVIVVVFAVCWLP LHVFNLL RDLDPHAIDPYAFGLVQLLCHWLMSSACYNPFYAWLHDSFREELRKLLVAW PRKIAPH GQNMTVSVVI
Source	in vitro E.coli expression system
Target Names	PRLHR
Protein Names	Recommended name: Prolactin-releasing peptide receptor Short name= PrRP receptor Short name= PrRPR Alternative name(s): G-protein coupled receptor 10 hGR3
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.