



Recombinant Rat Galanin receptor type 3 (Galr3)

Product Code	CSB-CF009225RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88626
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MADIQNISLDSPGSVGAVAVPVIFALIFLLGMVGNGLVLAVLLQPGPSAWQEPR STTDLF ILNLAVADLCFILCCVPFQAAIYTLDAWLFGAFVCKTVHLLIYLTMYASSFTLAAV SLDR YLAVRHPLRSRALRTPRNARAAGVGLVWLLAALFSAPYLSYYGTVRYGALELCV PAWEDAR RRALDVATFAAGYLLPVAVVSLAYGRTLCLFLWAAVGPAGAAAAEARRRATGR AGRAMLAV AALYALCWGPHHALILCFWYGRFAFSPATYACRLASHCLAYANSCLNPLVYSL ASRHFRA RFRRLWPCGRRRHRHHHRAHRALRRVQPASSGPAGYPGDARPRGWSMEP RGDALRGGGET RLTLSPRGPQ
Source	in vitro E.coli expression system
Target Names	Galr3
Protein Names	Recommended name: Galanin receptor type 3 Short name= GAL3-R Short name= GALR-3
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.