



Recombinant Rat Galanin receptor type 2 (Galr2)

Product Code	CSB-CF009224RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O08726
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MNGSGSQGAENTSQEGGSGGWQPEAVLVPLFFALIFLVGTVGNALVLAVLLR GGQAVSTT NLFILNLGVADLCFILCCVFPQATIYTLDDWVFGSLLCKAVHFLIFLTMHASSFTL AAVS LDRYLAI RYPLHSREL RTPRNALAAIGLIWGLALLFSGPYLSYYRQSQLANLTV C HPAWS APRRRAMDLCTFVFSYLLPVLVLSLTYARTLRYLWRTVDPVTAGSGSQRAKRK VTRMIII VAVLFCLCWMPHHALILCVWFGRFPLTRATYALRILSHLVSYANSCVNPIVYAL VSKHFR KGFRKICAGLLRPAPRRASGRVSI LAPGNHSGSMLEQESTDLTQVSEAAGPLV PPPALPN CTASSRTLDPAC
Source	in vitro E.coli expression system
Target Names	Galr2
Protein Names	Recommended name: Galanin receptor type 2 Short name= GAL2-R Short name= GALR-2
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.