



Recombinant Bovine Probable G-protein coupled receptor 173 (GPR173)

Product Code	CSB-CF704833BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5E9H8
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MANTTGEPEEVSGALSPPSAVAYVKLVLLGLIMCVSLAGNAILSLLVLKDRALH KAPYYF LLDLCLADGIRSAVCFPFVLASVRHGSSWTFALSCKIVAFMAVLFCFHAAFML FCISVT RYMAIAHHRFYAKRMTLWTCAAVICMAWTLVAMAFPPVFDVGTYKFIREEDQ CIFEHRY FKANDTLGFMLMLAVLMAATHAVYGKLLLFYRHRKMKPVQMMPAISQNWTF HGPGATGQ AAANWIAGFGRGMPPTLLGIRQNGHAASRLLGMDEVKGEKQLGRMFYAIT LLFLLLWS PYIVACYWRVFKACAVPHRYLATAVWMSFAQAAVNPIVCFLLNKDLKKCLRT HAPCWGT GGAPAPREPYCVM
Source	in vitro E.coli expression system
Target Names	GPR173
Protein Names	Recommended name: Probable G-protein coupled receptor 173
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.