



Recombinant Rat Probable G-protein coupled receptor 173 (Gpr173)

Product Code	CSB-CF888325RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JJH2
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MANTTGEPEEVSGALSLPSASAYVKLVLLGLIMCVSLAGNAILSLLVLKERALHK APYYF LLDLCLADGIRSAICFPFVLASVRHGSSWTFALSCKIVAFMAVLFCFHAAFMLF CISVT RYMAIAHHRFYAKRMTLWTCAAVICMAWTLVAMAFPPVFDVGTYKFIREEDQ CIFEHRY FKANDTLGFMLMLAVLMAATHAVYGKLLLFYRHRKMKPVQMMPAISQNWTF HGPGATGQ AAANWIAGFGRGMPPTLLGIRQNGHAASRLLGMDEVKGEKQLGRMFYAIT LLFLLLWS PYIVACYWRVFKACAVPHRYLATAVWMSFAQAAVNPIVCFLLNKDLKKCLRT HAPCWGT GGAPAPREPYCVM
Source	in vitro E.coli expression system
Target Names	Gpr173
Protein Names	Recommended name: Probable G-protein coupled receptor 173 Alternative name(s): Super conserved receptor expressed in brain 3
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.