



Recombinant Schizosaccharomyces pombe Protein ERD1 homolog 1 (SPAC227.01c, SPAPB21F2.04c)

Product Code	CSB-CF892227SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UTD8
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MALIEDLNHFINYFPLVLRLLFFLVVFGLYSFTLILHLLVINRVDVFSLLHTPLPVNR SQQ ANAPLWQLSFSL SILGTLLFVIAESLYLISGSDELAYVPVFIFGVIVFMPVHKFWF FQRK VFTRQCLRILGGSYRPDYKFPDVIFSDLLTSYSRVIADLWLAGAILIYVTDSPNN SHRKQ YENEVIMSMIAAYPYAIRFRQCLIERSSADNSSDKFWSTLNSIKYFTAFFPAIFLGI FAKK RFSFLWFLWNTSSAINSTYSFWWDVSMDSLPPFFKQPLSIQNWKFGVRRLLFP TTFFAVVS AIDFVLRMAWVVRVLPEHQSAFFTTDFGIFIMQFLEVFRRCVWVFFRIEAEASK SLAYVN ISDRSDIPTIHPD
Source	in vitro E.coli expression system
Target Names	erd1
Protein Names	Recommended name: Protein ERD1 homolog 1
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.