



Recombinant Mouse P2Y purinoceptor 2 (P2ry2)

Product Code	CSB-CF017334MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35383
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAADLEPWNSTINGTWEGDELGYKCRFNEDFKYVLLPVSYGVVVCVLGLCLNV VALYIFLC RLKTWNASTTYMFHLAVSDSLYAASLPLLVYYYARGDHWPFSTVLCKLVRFLF YTNLYCS ILFLTCSVHRCLGVLRPLHSLRWGRARYARRVAAVVWVLVLACQAPVLYFVTT SVRGTR ITCHDTSARELFSHFVAYSSVMLGLLFAVPFSVILVCYVLMARRLLKPAYGTTG GLPRAK RKSVRTIALVLAVFALCFLPFHVTRTLYYSFRSLDLSCHTLNAINMAYKITRPLAS ANSC LDPVLYFLAGQRLVRFARDAKPTEPTPSPQARRKLGLHRPNRTVRKDLSVSS DDSR RTE STPAGSETKDIRL
Source	in vitro E.coli expression system
Target Names	P2ry2
Protein Names	Recommended name: P2Y purinoceptor 2 Short name= P2Y2 Alternative name(s): ATP receptor P2U purinoceptor 1 Short name= P2U1 Purinergic receptor
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.