



Recombinant Guinea pig P2Y purinoceptor 1 (P2RY1)

Product Code	CSB-CF017326GU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59902
Product Type	Transmembrane Protein
Immunogen Species	Cavia porcellus (Guinea pig)
Sequence	MTEVLWPAAPNGTDA AFLASPGFHWGNSTATSTAAAAAPFRCALTKTGFQFY YLPVYIV VFIIGFLGNSIAIWMFVFHMKPWSGISVYMFNLALADFLYVLTLPALIFYFYNKTN WIFG DAMCKLQRFIFHVNLVYGSILFLTCISAHRYSGVVYPLKSLGRLKKKNAVYISVLV WLIVV VAISPILFYSGTGIRKNKTITCYDTSDEYLRSYFIYSMCTTVAMFCVPLVLILGC YGLI VRALIYKDLDNSPLRRKSIYLVIIIVLTVFAVSYIPFHVMTMNLRARLDFQTPEMC TFND RVYATYQVTRGLASLNSCVDPILYFLAGDTRRRRLSRATRKASRRSEANLQSK SEDMTLN ILSEFKQNGDTSL
Source	in vitro E.coli expression system
Target Names	P2RY1
Protein Names	Recommended name: P2Y purinoceptor 1 Short name= P2Y1 Alternative name(s): ATP receptor Purinergic receptor
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.