



# Recombinant Schizosaccharomyces pombe Uncharacterized transporter C22E12.01 (SPAC22E12.01, SPAC890.09)

<b>Product Code</b>	CSB-CF611958SXV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q10354
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Sequence</b>	METDPILKVEREPEFTPSPAPEAIVKIDIEQGVAPVSQNEPDRPQHWITRVVIIVLIVLAW YFFSLLLMMNKWIFSESKMDFQFPLFLSSCQMLVQMGMFAKLILAFPRYQPN KKDNFSW LEYFYRAGICALVTGLDIGLSNASLETITLSFYTMCRSSILIFVFFFSVIFRIEMFD WIL LCITLVISAGVVLMTATETQFVLSGFLVMASVLSGLRWALTQKLLLDHPWTS NPFTSL FALTPLMFLFLLVAGLIFEGPVRFIESPAWKEFGPFMSVILVPGTLAFFMVASE FGLIQ KTSIVTLVCGILKEIITIIASTLFYHDILLPINIVGLVITLCGIGVYNYRITKGNKKE AEKEVEYIVLNENA
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SPAC22E12.01
<b>Protein Names</b>	Recommended name: Uncharacterized transporter C22E12.01
<b>Expression Region</b>	1-374
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.