



Recombinant Mouse 5-hydroxytryptamine receptor 1D (Htr1d)

Product Code	CSB-CF726741MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61224
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MSPPNQSLEGLPQEASNRSLNVTGAWDPEVLQALRISLVVLSVITLATVLSNA FVLTTI LLTKKLHTPANYLIGSLATDLLVSILVMPISIAYTTRTWNFGQILCDIWWSSDIT CCT ASILHLCVIALDRYWAITDALEYSKRRTAGHAAAMIAAVWIISICISIPPLFWRQAT AHE EMSDCLVNTSQISYTIYSTCGAFYIPSILLIILYGRIYVAARSRLNPPSLYGKRFT TAQ LITGSAGSSSLCSLNPSLHESHTHTVGSPLFFNQVKIKLADSILERKRISAARERK ATKTL GIILGAFIICWLPFFVVSLVLPICRDSCWIHPALFDFFTWLGYLNSLINPVIYTVFN EDF RQAFQKVVFHFRKIS</p>
Source	in vitro E.coli expression system
Target Names	Htr1d
Protein Names	Recommended name: 5-hydroxytryptamine receptor 1D Short name= 5-HT-1D Short name= 5-HT1D Alternative name(s): Serotonin receptor 1D
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.