



Recombinant Serpentine receptor class alpha-39 (sra-39)

Product Code	CSB-CF527677CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O62420
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MEGCETRGQLYYHPAYLFAIIFQGTSGYLTFPILFIILKKYAFKVYYHPNLVFLMM LNVV SCLLLGGLTAWSATNFFLNLMIKPSPCDLLLKTYFCISKIRGTFLFAFCLVTTSHA GILLE RSWATYSVKNYERQGRALGTILAIVAVAVAATAIFILLLPEDGEEMITTCLTFSAS KSIG SRVYVMFFVQLLLLDAVISIVHLVLYRYNKNIDKHGSTLSSEQFQRNENVKTLKQ VTPLLI LSNVTIGVYIFIVSVFRLYKDYLPNWEIIAANLFIMPHMPFMFTSLILVELWLG KKRQ DRIHKDMMASQDVPLDDQFHAMENWDVHFQSRLKAQARANRKPLSANWTS STLITKFKK RQTTVAVISIPSPS
Source	in vitro E.coli expression system
Target Names	sra-39
Protein Names	Recommended name: Serpentine receptor class alpha-39 Short name= Protein sra-39
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.