



# Recombinant Human G-protein coupled receptor family C group 5 member B (GPRC5B)

<b>Product Code</b>	CSB-CF868367HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9NZH0
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	ENASTSRGCGLDLLPQYVSLCDLDAIWGIVVEAVAGAGALITLLLMLILLVRLPFI KEKE KKSPVGLHFLFLLGLTGLFGLTFAFIIQEDETICSVRRFLWGVLFALCFSCLLSQ AWRVR RLVRHGTGPAGWQLVGLALCLMLVQVIIAVEWLVLTVLRDTRPACAYEPMDFV MALIYDM VLLVVTLGLALFTLCGKFKRWKLNQAFLLITAFLSVLIWVAWMTMYLFGNVKLQ QGDAWN DPTLAITLAASGWVVFVIFHAPEIHCTLLPALQENTPNYFDTSQPRMRETAFEED VQLPR AYMENKAFSMDEHNAALRTAGFPNGSLGKRPSGSLGKRPSAPFRSNVYQPT EMAVVLNGG TIPTAPPSHTGRHLW
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	GPRC5B
<b>Protein Names</b>	Recommended name: G-protein coupled receptor family C group 5 member B Alternative name(s): A-69G12.1 Retinoic acid-induced gene 2 protein Short name= RAIG-2
<b>Expression Region</b>	29-403
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.