



Recombinant Human Probable G-protein coupled receptor 27 (GPR27)

Product Code	CSB-CF885717HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NS67
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MANASEPGGSGGGEEAALGLKLATLSLLLCVSLAGNVLFALLIVRERSLHRAPY YLLLDL CLADGLRALACLPVMLAARRAAAAAGAPPGALGCKLLAFLAALFCFHAAFLLL GQGVTR YLAIHHRFYAERLAGWPCAAMLVCAAWALALAAFPVLDGGGDDDEDAPCA LEQRPDGA PGALGFLLLAVVVGATHLVYLRLFFIHDRRKMRPARLVPAVSHDWTFFHGPG ATGQAAA NWTAGFGRGPTPPALVGIRPAGPGRGARRLLVLEEFKTEKRLCKMFYAVTLLF LLWGPY VVASYLRVLRPGAVPQAYLTASVWLTFQAQAGINPVVCFLFNREL RDCFRAQF PCCQSPR TTQATHPCDLKGIGL
Source	in vitro E.coli expression system
Target Names	GPR27
Protein Names	Recommended name: Probable G-protein coupled receptor 27 Alternative name(s): Super conserved receptor expressed in brain 1
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.