



Recombinant Human Neuropeptide Y receptor type 4 (PPYR1)

Product Code	CSB-CF018593HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50391
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNTSHLLALLLPKSPQGENRSKPLGTPYNFSEHCQDSVDVMVFIVTSYSIETVV GVLGNL CLMCVTVRQKEKANVTNLLIANLAFSDFLMCLLCQPLTAVYTIMDYWIFGETLC KMSAFI QCMSVTVSILSLVLVALERHQLIINPTGWKPSISQAYLGIVLIWVIACVLSLPFLA NSIL ENVFHKNHSKALEFLADKVVCTESWPLAHHRTIYTTFLLLFQYCLPLGFILVCYA RIYRR LQRQGRVFKGTYSRLRAGHMKQVNVVLVVMVVAFAVLWLPLHVFNSLEDWH HEAIPICHG NLIFLVCHLLAMASTCVNPFYIGFLNTNFKKEIKALVLTCCQSAPLEESEHLPLS TVHTE VSKGSLRLSGRSNPI
Source	in vitro E.coli expression system
Target Names	NPY4R
Protein Names	Recommended name: Neuropeptide Y receptor type 4 Short name= NPY4-R Alternative name(s): Pancreatic polypeptide receptor 1 Short name= PP1
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.