



Recombinant Inner membrane transport permease ybhS (ybhS)

Product Code	CSB-CF360304EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AFQ4
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MSNPILSWRRVRALCVKETRQIVRDPSSWLIAVVIPLLLLFIFGYGINLDSSKLRV GILL EQRSEAALDFHTMTGSPYIDATISDNRQELIAKMQAGKIRGLVVIPVDFAEQM ERANAT APIQVITDGSEPNTANFVQGYVEGIWQIWQMQRAREDNGQTFEPLIDVQTRYWF NPAAISQ HFIIPGAVTIIMTVIGAILTSLVVAREWERGTMEALLSTEITRTELLLCKLIPYYFLG ML AMLLCMLVSVFILGVPYRGSLILFFISSLFLLSTLGMGLLISTITRNQFNAAQVAL NAA FLPSIMLSGFIFQIDSMPIRAVITYIIPARYFVSTLQSLFLAGNIPVVLVNVLFLLI AS AVMFIGLTWLKTKRRLD
Source	in vitro E.coli expression system
Target Names	ybhS
Protein Names	Recommended name: Inner membrane transport permease ybhS
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.