



Recombinant *Macaca mulatta* Apelin receptor (APLNR)

Product Code	CSB-CF001910MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O97666
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca mulatta</i> (Rhesus macaque)
Sequence	<p>MEEGGDFDNYYGADNQSECEYTDWKSSGALIPAIYMLVFLLGTTGNGLVLWT VFRSSREK RRSADIFIASLAVADLTFVVTLPPLWATYTYRDYDWPFGTFSCKLSSYLIFVNM YASVFCL TGLSFDRLAIVRPVANARLRLRVSGAVATAVLWVLAALLAMPVMVFRRTTGDL ENTTKVQ CYMDYSMVATVSSDWAVEVGLGVSSTTVGFVVPFTIMLTCYFFIAQTIAGHFR KERIEGL RKRRRLLSIIVVLVVTFALCWMPYHLVKTLYMLGSLHWP CDFDLFLMNVPYC TCISYV NSCLNPFLYAFFDPRFRQACTSMLCCGQSRCAGTSHSSSGEKSASYSSGHS QGPGPNMGK GGEQMHEKSIPYSQETLVVD</p>
Source	in vitro E.coli expression system
Target Names	APLNR
Protein Names	Recommended name: Apelin receptor Alternative name(s): Angiotensin receptor-like 1 G-protein coupled receptor APJ
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.