



Recombinant Drosophila melanogaster Putative gustatory receptor 39a, isoform D (Gr39a)

Product Code	CSB-CF348430DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58958
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MKRNAFEELRVQLRTLKWLGVLRFTIDFNKCLVRENASEERSAWLYLIGVVGIT CSLIVY STYFPSHFIMGKHNTTGNCYALINIRSCSIVTMLIYTQLYIQRFRFVALLQSILRF NQIS GSHREEGRFAFYYYTHLSLLIICMLNYAYGYWTAGVRLTTIPIYLLQYGFSYLFL GQVVV LFACIQQILLSILKYYNQVVLKNIKSSKESREFYYNFCKYNQVIWLSYTEINHCFG LLLL LVTGLILLITPSGPFYLVSTIFEGRFRQNWQFSLMSFTAILWSLPWIVLLVLAMG RNDVQ KEANKTAKMLTKVPRTGTGLDRMIEKFLKLNLRQKPILTAYGFFALDKSTLFLKF TAIFT YMVILVQFKEMENSTKSINKF
Source	in vitro E.coli expression system
Target Names	Gr39a
Protein Names	Recommended name: Putative gustatory receptor 39a, isoform D
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.