



Recombinant *Drosophila melanogaster* Putative gustatory receptor 39a, isoform C (Gr39a)

Product Code	CSB-CF349275DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58957
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MDFQPGELCAYYRLCRYLGIFCIDYNPTKKKFRLRRSVLCYIVHFALQAYLVGCI SVMVT YWRRCFKSELTTTGNHFDRLVMVIALGILVVQNAWLIWLQAPHLRIVRQIEFYR RNHLAN VRLLLPKRLLWLIATNVVYMANFIKTCIFEWLTASRLFVITSLGFPLRYLVTSFT MGT YFCMVHIVRLVLDWNQSQINAIIDESADLKMTSPNRLRLRVCLEMHDRLMLLCN DEISLV YGFIAWLSWMFASLDVTGVYILTMVIQTKKSIVLKLITNVVWLSPTFMTCAASFM SNRVT IQANKTAKMLTKVPRTGTGLDRMIEKFLLNLRQKPILTAYGFFALDKSTLFLKF TAIFT YMVILVQFKEMENSTKSINKF
Source	in vitro <i>E.coli</i> expression system
Target Names	Gr39a
Protein Names	Recommended name: Putative gustatory receptor 39a, isoform C
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.