



Recombinant Gorilla gorilla gorilla Probable G-protein coupled receptor 34 (GPR34)

Product Code	CSB-CF747689GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6XCF2
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	MRSHTITMTTTSVSSWPYSSHRMRFITNHSDQPPQNFSA TPNVTTCPMDEKLL STVLTTST YSVIFIVGLVGNIIALYVFLGIHRKRNSIQIYLLNVAIADLLLIFCLPFRIMYHINQNK W TLGVILCKVVGTLFYMNMYISIILLGFISLDRIYIKINRSIQQRKAITTKQSIYVCCIV WM LALGGFLTMILTLLKGGHNSTMCFHYRDKHNAKGEAIFNFILVVMFWLIFLLIILS YIK IGKNLLRISKRRSKFPNSGKYATTARN SFIVLIIFTICFVPYHAFRFIYISSQLNVSS CY WKEIVHKTNEIMLVLSSFNSCLDPV MYFLMSSNIRKIMCQLLFRRFQGEPSRSE STSEFK PGYSLHDT SVAVKIQSSSKST
Source	in vitro E.coli expression system
Target Names	GPR34
Protein Names	Recommended name: Probable G-protein coupled receptor 34
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.