



Recombinant Dictyostelium discoideum Putative UDP-sugar transporter DDB_G0278631 (DDB_G0278631)

Product Code	CSB-CF681792DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q54YK1
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	MKLPNDIYIPLNSIQFKIHNDINIKYDTEEDNNNNNNNNNNNNRYIDDNEKMKI SEDES IFKSINNDNKNQQNDNIINSFNRSNSFLFNIINFYKSPRFLSGGAVLLYVMTTY DFSAS NFLLFNQMVVTIVILHILKHFNILKINTNFEKETIKKLMPLSFCYIINVLLGLDSLKQ LN IPMYSALKRLVAVVILVMEYFILKKVSPPKIIASVVVMVIGAVVAGITDLSFNLSLGY SLV LLSCIFQASYLIYVKKVASNMSTYDMLYNSVLSLPITIFLMIVNQEIEYFQTFEHL YDS SFQAYFILSIFLGFFLNFCIFFCTSVNSPLTTSVTGQVKNIASTIIGAMVFNDIIHPI N ILGLIINIIGSIWYSFLKLTSK
Source	in vitro E.coli expression system
Target Names	DDB_G0278631
Protein Names	Recommended name: Putative UDP-sugar transporter DDB_G0278631
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.