



Recombinant Drosophila melanogaster Putative odorant receptor 94b (Or94b)

Product Code	CSB-CF893733DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VCS8
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MESTNRLSAIQTLLVIQRWIGLLKWENEGEDGVLTWLKRIYPFVLHLPLTFTYIA LMWYE AITSSDFEEAGQVLYMSITELALVTKLLNIWYRRHEAASLIHELQHDPANLNRNS EEIKF WQQNQQRNFKRIFYWYIWGSLFVAVMGYISVFFQEDYELPFGYYVPFEWRTRE RYFYAWGY NVVAMTLCCLSNILLDTLGCYFMFHASLFRLLGMRLEALKNAAEKARPELRRRI FQLHT KVRRLTRECEVLVSPYVLSQVVFSAFIICFSAYRLVHMGFKQRPGLFVTTVQFV AVMIVQ IFLPCYYGNELTFHANALTNSVFGTNWLEYSVGTRKLLNCYMEFLKRPVKVRA GVFFEIG LPIFVKTINNAYSFFALLLKISK
Source	in vitro E.coli expression system
Target Names	Or94b
Protein Names	Recommended name: Putative odorant receptor 94b
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.