



Recombinant Organic solute transporter alpha-like protein C01B12.4 (C01B12.4)

Product Code	CSB-CF517701CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O17204
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MEIVKTIIPHNRSYIEPPISATEWLANMSVMHVSCLTIACVFVAITFLSSFFHLFF VLK YVSNERIRNDMYALIFMFPITTFASLVGMFIPRAAIFLYAVSLVYFMFTLFIMVTLL FNI FGGRQEMSAYLLQRNIRVNFTVPPLCFFKFLPTVESTDQNLRRIEWLVFQTPII RTLLEL VSVVVSMEQEGRRRESVWFVFSQLMALLSMCIAFYGCYVMVPLGREKHAPYRF DFLFRTCD IAQCIYTIQKFVFEFAAAVGLITSDRYLPAAAKALWWASFMCTWEMMLLSALCS YCLRPA KCKFFDLYPGNDMPALSARDGSNSRVPSFSRRLSIEYEPRIAGVMLEPPSRSS LSITPRD KIEDPTTVSYFADNFDSLSQIQGQ
Source	in vitro E.coli expression system
Target Names	osta-1
Protein Names	Recommended name: Organic solute transporter alpha-like protein C01B12.4
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.