



Recombinant Bovine Prokineticin receptor 2 (PROKR2)

Product Code	CSB-CF018749BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8SPN1
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MAAQNGNASFPANFSIPQEHASSLPFNFSYDDYDLPLDEDEDMTKTQTFFAAK IVIGVAL VGIMLTCGIGNFVFITALTRYKKLRNLNLLIANLAISDFLVAIICCPFEMDYYVVH QLS WEHGHVLCACINYLRVSLYVSTNALLAIAIDRYLAIVHPLKPRMNYQTASFLIAL VWMV SILISIPSAYFTKETVLFIVKNQKKIFCGQVWPVDQQLYYKSYFLFVFGIEFLGPV FTMT LCYARISRELWFKAVPGFQTEQIRKRLRCRRKTVLVLMLCILTAYVLCWAPFYGF TIVRDF FPTVFVKEKHLYLTAFYVVECIAMSNMINTVCFVTVKNSTMKYFKKMLLLHWR PSHHGSK SSADLDLKTSRLPATEEVDCIRLK
Source	in vitro E.coli expression system
Target Names	PROKR2
Protein Names	Recommended name: Prokineticin receptor 2 Short name= PK-R2 Alternative name(s): G-protein coupled receptor 73-like 1 G-protein coupled receptor I5E
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.