



Recombinant Mouse Probable G-protein coupled receptor 88 (Gpr88)

Product Code	CSB-CF887640MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9EPB7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MTNSSSTSTSTTTGGSLLLLCEEEESWAGRRIPVSLLYSGLAIGGTLANGMVIY LVSSFR KLQTTSNAFIVNGCAADLSVCALWMPQEAVLGLLPSGSAEPPGDWDGGGGS YRLLRGGLL GLGLTVSLLSHCLVALNRYLLITRAPATYQVLYQRRHTVGMLALSWALALGLVL LLPPWA PKPGAEPQVHYPALLAAGALLAQTALLLHCYLGIVRRVRVSVKRVSVLNFHLL HQLPGC AAAAAFPAAPHAPGPGGAAHPAQPLPAALQPRRAQRRLSGLSVLLLCCV FLLATQPL VWVSLASGFSLPVPWGVQAASWLLCCALSALNPLLYTWRNEEFRRSVRSVLP GVGDAAAA AAAATAVPAMSQAQLGTRAAGQHW
Source	in vitro E.coli expression system
Target Names	Gpr88
Protein Names	Recommended name: Probable G-protein coupled receptor 88 Alternative name(s): Striatum-specific G-protein coupled receptor
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.