



Recombinant *Drosophila melanogaster* Putative odorant receptor 98b (Or98b)

Product Code	CSB-CF895586DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VAW0
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MLTDKFLRLQSALFRLGLLELLHEQDVGHRYPWRSICCILSVASFMPPLTIAFGL QNVQNV EQLTDSLCSVLVDLLALCKIGLFLWLYKDFKFLIGQFYCVLQTETHHTAVAEMIVT RESRR DQFISAMYAYCFITAGLSACLMSPLSMLISYQRTGELQPKFPFSPVYPWDNMK LSNYIIS YFWNVCAALGVALPTVCVDTLFCSLSHNLCALFQIARHKMMHFEGRNTKETHE NLKHVFQ LYALCLNLGHFLNEYFRPLICQFVAASLHLCVLCYQLSANILQPALLFYAAFTAA VVGQV SIYCFCGSSIHSECQLFGQAIYESSWPHLLQENLQLVSSLKIAMMRSSLGCPID GYFFEA NRETLITVSKAFIKVSKKTPQVND
Source	in vitro <i>E.coli</i> expression system
Target Names	Or98b
Protein Names	Recommended name: Putative odorant receptor 98b
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.