



Recombinant Bovine Neuropeptide Y receptor type 2 (NPY2R)

Product Code	CSB-CF016036BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79113
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MKMGPLGAEADENQTV EEMKVDQFGPGHTTLPGELAPDSEPELIDSTKLIEVQ VVLILAY CSIILLGVIGNSLVIHVVIKFKSMRTVTNFFIANLAVADLLVNTLCLPFTLTYTLMG EWK MGPVLCHLVPYAQGLAVQVSTITLTVIALDRHRCIVYHLESKISKQISFLIIGLAW GVSA LLASPLAIFREYSLIEIIPDFEIVACTEKWPGEEKGIYGTIYSLSSLLILYVLPGLIIS F SYTRIWSKLNHVSPGAAHDHYHQRRQKTTKMLVCVVVFAVSWLPLHAFQL AVDIDSHV LDLKEYKLIFTVFHIIAMCSTFANPLLYGWMNSNYRKAFLSAFRCEQRLDAIHSE VSVTF KAKKHLQVTKNNGPNDSTETTIV
Source	in vitro E.coli expression system
Target Names	NPY2R
Protein Names	Recommended name: Neuropeptide Y receptor type 2 Short name= NPY2-R Alternative name(s): NPY-Y2 receptor Short name= Y2 receptor
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.