



Recombinant Papio hamadryas Proteinase-activated receptor 1 (F2R)

| | |
|--------------------------|--|
| Product Code | CSB-CF007924ERJ |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P56488 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Papio hamadryas (Hamadryas baboon) |
| Sequence | SFLLRNPNDKYEPFWEDEEKNESGLTEYRLVSINKSSPLQKPLPAFISEDASGY LTSSWL TLFVPSVYTGVFVVS L P V N I M A I V V F I L K M K V K P A V V Y M L H L A T A D V L F V S V L P F K I S Y Y L S G S D W Q F G S E L C R F V T A A F Y C N M Y A S I L L M T V I S I D R F L A V V Y P M Q S L S W R T L G R A S F T C L A I W A L A I A G V V P L L L K E Q T I Q V P G L N I T T C H D V L N E T L L E G Y Y A Y Y F S A F S A V F F F V P L I I S T V C Y V S I I R C L S S S T V A N R S K K S R A L F L S A A V F C I F I I C F G P T N I L L I A H Y S F L S H T S T T E A A Y F A Y L L C V C V S S I S C C I D P L I Y Y Y A S S E C Q R Y V Y S I L C C K E S S D P S S S N S S G Q L M A S K M D T C S S N L N N S I Y K K L L T |
| Source | in vitro E.coli expression system |
| Target Names | F2R |
| Protein Names | Recommended name: Proteinase-activated receptor 1 Short name= PAR-1 Alternative name(s): Thrombin receptor |
| Expression Region | 42-425 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |