



Recombinant Mouse Neuropeptide Y receptor type 2 (Npy2r)

Product Code	CSB-CF016036MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P97295
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MVLKMGVPGAEADENQTVEVKVEPYGPGHTTPRGELPPDPEPELIDSTKLVE VQVILILA YCSIILLGVVGNLSLVHVVVIFKFSMRTVTNFFIANLAVADLLVNTLCLPFTLTITLM GEW KMGPVLCHLVPYAQGLAVQVSTITLTVIALDRHRCIVYHLESKISKRISFLIIGLA WGIS ALLASPLAIFREYSLIEIIPDFEIVACTEKWPGEKSVYGTVYSLSTLLILYVPLPLGI IS FSYTRIWSKLRNHVSPGAASDHYHQRRHKMTKMLVCVVVFAVSWLPLHAFQ LAVDIDSH VLDLKEYKLIFTVFHIIAMCSTFANPLLYGWMNSNYRKAFLSAFRCEQRLDAIHS EVSMT FKAKKNLEVKKNNGPTDSFSEATNV
Source	in vitro E.coli expression system
Target Names	Npy2r
Protein Names	Recommended name: Neuropeptide Y receptor type 2 Short name= NPY2-R Alternative name(s): NPY-Y2 receptor Short name= Y2 receptor
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.