



Recombinant Drosophila melanogaster Putative odorant receptor 82a (Or82a)

Product Code	CSB-CF307306DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P82986
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MGRLFQLQEYCLRAMGHKDDMDSTDSTALSLKHISLIFVISAQYPLISYVAYN RNDMEK VTACLSVVFTNMLTVIKISTFLANRKDFWEMIHFRFRKMHEQSASHIPRYREGLD YVAEAN KLASFLGRAYCVSCGLTGLYFMLGPIVKIGVCRWHGTTCDKELPMPMKFPFND LESPGYE VCFLYTVLVTVVVVAYASAVDGLFISFAINLRAHFQTLQRQIENWEFPSSEPDT QIRLKS IVEYHVLLLSSLRKLRSIYTPVMGQFVITSLQVGVIIYQLVTNMDSVMDLLLYAS FFGS IMLQLFIYCYGGEIKAESLQVDTAVRLSNWHLASPKTRTSLSLIILQSQKEVLIRA GFF VASLANFVGICRTALSLITLIKSIE
Source	in vitro E.coli expression system
Target Names	Or82a
Protein Names	Recommended name: Putative odorant receptor 82a
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.