



Recombinant Human Adiponectin receptor protein 2 (ADIPOR2)

Product Code	CSB-CF768214HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q86V24
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNEPTENRLGCSRTPEPDIRLRKGHQLDGTRRGDNDSHQGDLEPILEASVLS SHHKKSSE EHEYSDEAPQEDEGFMGMSPLLQAHAMEKMEEFVCKVWEGRWRVIPHDVL PDWLKDNDF LLHGHRPPMPSFRACFKSIFRIHTETGNIWTHLLGCVFFLCLGIFYMFRPNISFV APLQE KVVFGFLFFLGAILCLSFSWLFHTVYCHSEGVSRLFSKLDYSGIALLIMGSFVPWL YYSFY CNPQPCFIYLIVICVLGIAAIVSQWDMFATPQYRGVVRAGVFLGLGLSGIIPTLHY VISE GFLKAATIGQIGWLMLMASLYITGAALYAARIPERFFPGKCDIWFHSHQLFHIFV VAGAF VHFHGVSNLQEFRFMIGGGCSEEDAL
Source	in vitro E.coli expression system
Target Names	ADIPOR2
Protein Names	Recommended name: Adiponectin receptor protein 2 Alternative name(s): Progestin and adipoQ receptor family member II
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.